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Combining Flexibility and Ultrahigh Performance in Mass Spectrometry

Michael L. Easterling Ph.D.
Bruker Daltonics, Inc. Billerica, MA

As a result of the ongoing expansion in the types of mass spectrometric tools and workflows available to the scientific community, one of the toughest challenges which faces the modern user is to properly outfit themselves with the systems needed to efficiently answer analytical questions. From selection of an appropriate sample cleanup mode to an ionization strategy, to a detection technique, and a variety of post processing options, the parameter space for technique selection is seemingly infinite. The power derived from this flexibility has created a unique opportunity for expansion in mass spectrometry to a variety of application areas that have not been conventionally considered.

Historically, Fourier transform ICR mass spectrometry has been viewed as an “event horizon” platform for the early adaptation or development of many tools used today in conventional mass spectrometry, to include; electron induced dissociation techniques, electrospray ionization, top-down MS, and HDX. Many of these tools developed for the life sciences either had their genesis in, or were refined on FTMS platforms. This tradition is carried on today with the incorporation of MALDI imaging and other advanced or emerging techniques into FTMS instrumentation where these technologies can utilize the unique capabilities and spectral figures-of-merit found in FTMS to answer questions in a fashion that complements other MS strategies. While many types of mass spectrometers can provide access to workflows attributed to FTMS, none can leverage the additional dimension of ultrahigh performance that makes FTMS the most powerful and flexible tool in the mass spec arsenal.

This talk will examine how ultrahigh performance can provide definitive molecular formula determination and identification for a variety of application focus areas. Specifically, we will look at high resolving power (>250,000 FWHM) mass spectrometry applied to a variety of workflows. In each case, the “extra” chemical information encoded in the isotopic peaks serves as a unique signature to the atomic composition of the target molecule. While methodologies based on analysis of the geometric features of unit-resolved isotopic distributions have become fashionable to overcome the degeneracy associated with mass based elemental assignment, they cannot definitively exclude the possibility of alternative atomic combinations or mixtures. Measuring the individual atomic contribution to each of these peaks can uniquely provide complete confidence formula elucidation, even for complex mass spectra

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Challenges in Metabolomics and metabolism studies addressed by UHR-TOF-MS technology combined with unique software capabilities

Aiko Barsch Ph.D.,
Bruker Daltonics, Bremen, Germany

An MS system's performance in terms of mass accuracy, resolution, dynamic range, sensitivity and MS/MS performance are often pushed to their limits with complex samples originating from Metabolomics or drug Metabolism studies. The micrOTOF and maXis instrument series provide an unparalleled combination of the necessary MS performance features required to analyze these highly complex samples. When utilized in conjunction with tailored software solutions, this potent combination perfectly addresses the needs in Metabolomics and MetaboliteID for accurate and comprehensive small molecule analysis.

One crucial advantage of ultra high resolution UHR-TOF technology is the speed of analysis by coupling U-HPLC systems without compromising on the performance in terms of mass accuracy or resolution even at high acquisition rates. Overall, the sample throughput is increased addressing the need for handling of large sample numbers.

Presently, identification of unknown metabolites is often defined as the major bottleneck in Metabolomics and likewise represents a daily challenge in drug metabolism studies. With increasing molecular mass of a compound, the number of possible molecular formulae increases exponentially. Even a mass accuracy of 0.1 ppm is not sufficient for an unambiguous formula identification for m/z values above 500. The unique SmartFormula 3D software tool extends formula generation to higher m/z values by combining mass accuracy, isotopic pattern, adduct and neutral loss information from MS and MS/MS data. An interactive link between SmartFormula 3D results, mass spectra and molecular structures is provided by the Fragment Explorer. This integrated tool is based on well-known ChemDraw™ technology for chemical structures, and has been especially designed for faster interpretation of MS/MS data.

Several examples will be presented which highlight the different aspects and benefits of applying UHR-MS technology for metabolomics and metabolism studies.

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Hyphenation of NMR, MS and Chromatography.
Manfred Spraul Ph.D.,
Bruker Biospin, Inc. Germany

NMR over the last years has strongly moved into mixture analysis, this can be either directly by looking into biofluids, food materials or extracts of any kind. Such NMR has developed into one of the main 2 technologies for metabolomics analysis, the other being mass spectrometry.

The strength of NMR is its unmatched reproducibility and transferability and its high dynamic range under full quantification. On the other hand MS has very high sensitivity and advantages when looking into large molecules in mixtures. Therefore the 2 technologies integrated for direct mixture analysis produce a lot of synergy. Therefore an integrated system has been developed, called the Metabolic Profiler, this system also contains a liquid handler for sample preparation and transfer. The system allows to measure NMR and LC-MS (preferably ultrahigh pressure LC) in a synchronized mode, this guarantees high quality results on integrated statistics like the heterostatistical correlation spectroscopy. Examples are given for biofluid and food analysis.

Another integrated system, which can be built up by expanding the Metabolic Profiler is the LC-SPE-NMR/MS system. In this case the task is to isolate compounds from a mixture to do structure confirmation or elucidation using NMR and MS. The system contains an LC-NMR interface, preferably a post column solid phase extraction unit. The system operates in such a way, that the MS (or MS + UV) information is used to decide on which peaks are interesting for post-column trapping and subsequent NMR analysis. With this technology NMR sensitivity can be boosted substantially, especially when transferring the trapped content into a small volume cryogenic probe. In addition multiple trapping allows to further enhance the NMR sensitivity. All this is possible under full automation. The system is introduced and examples are given.

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TOP-DOWN PROTEOMICS, A MULTI-PLATFORM PRIMER

Michael L. Easterling Ph.D.
Bruker Daltonics, Inc. Billerica, MA

Although the idea of “top-down” processing intact proteins has been around for a some time, only recently have software and hardware advances allowed this thesis to be pursued as a viable workflow. While there are many differing opinions about the most efficient way to gather information from top-down protein work, the most conventional misconception is that it requires the use of high performance instrumentation. Indeed, almost any kind of mass spectrometer capable of interfacing either MALDI or ESI can be used to generate data from intact proteins.

The challenge for the modern mass spectrometrist is to determine what kind of “hammer” to apply to the problem at hand. “Lower end” analyzers such as ion traps that offer lower resolving power and mass accuracy have been given new life by recent advances in the combination of advanced ion fragmentation techniques such as electron transfer dissociation (ETD) and increased analytical performance, and are unquestionably capable of performing top-down quality control (QC) roles.

Time-of-flight instrumentation coupled with electrospray ionization (ESI) has recently seen architectural improvements that allow resolving powers in excess of 50,000 allowing charge state and, ergo, mass determination of larger proteins even under fast chromatographic conditions. MALDI-TOF has seen recent methodological advances that allow protein terminal sequencing that rivals standard Edmond sequencing techniques in terms of sequence coverage.

Finally, the MS arsenal is completed with Fourier Transform (FTMS) based instruments which remain the most flexible and evaluative form of analyzer for the top-down protein chemists. This discussion will focus on how each of these platforms play an integral role in advancing top-down strategies and how each can be used in an integrated whole protein strategy to complement bottom-up workflows.

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